

neuropad[®]

Diagnostic Test for Sudomotoric Dysfunction (Anhydrosis)
and Test for the early detection of the Diabetic Foot Syndrome



neuropad[®] is a diagnostic tool to test sudomotoric functions by way of a colour changing pad which is applied at the sole of the foot (MTK I / II). This innovative test makes an objective, uncomplicated and economic examination of the diabetic patient possible and offers new therapeutic options without being invasive. Timely detection and intervention strategies can slow or even reverse the progress of some forms of neuropathy.

The test does not require the support of the patient and detects a neuropathy before the symptoms and signs of sensory and thermal disturbances are present in most cases.

Intended to objectively document the presences or absence of sweat as an adjunction to clinical assessment of diabetic patients as well as a major warning signal for autoimmune neuropathy.

For further clinical data please see our website:

http://trigocare.com/mediapool/94/944453/data/neuropad_Study_Total_110902-1_srx.pdf

¹ Retail Pack contains 1 Alufoil with 2 Testpatches + 1 Instruction for use

² Retail 10 Pack contains 10 Single Retail Packs

³ Clinic Pack 20pcs. Box contains 20 Alufoils + 1 Instruction for use

Order Information:

Packaging	Item-N°	DU
Retail Single ¹ Pack	16900	1 pcs.
Retail Single ² 10 Pack	16910	10 pcs.
Clinic Pack ³ 20pcs. Box	16920	20 pcs.

Composition of Goods:

Backing material: foil, transparent, flexible and water resistant

Adhesive: acrylate, hypoallergenic

Compress:

Backing material: cellulose „Paraprint OL 60“, 100% viscose (with acrylate adhesive)

Impregnating agent: CoCl₂ solution 6,25% (+/-0,25%)

Peel off strips: Paper, white, one side siliconized

widths: 2 x (40mm x 39mm)

Pouche: OPP-foil (single strips), transparent

widths: 54mm / lengths: 99mm

Storage Instructions:

First in - First out. Store away from direct sunlight at ambient temperature and humidity.

Do not store next to heating and/or entrance doors.

Stretchable ca. (in %): inapplicable

Sterilization: none-sterile

Washable: inapplicable

CE : 

Custom-Tariff-N°: 3005.1000

Shelf Life: 3 years from date of manufacturing

Classification: I

Annex: I + VII

Rule: 1

Product for single use yes

TRIGOCare International GmbH
In den Wiesen 3
D-51674 Wiehl
Germany

Fon: +49 (0)22 62 – 72701-0
Fax: +49 (0)22 62 – 7248-48
www.trigocare.com
info@trigocare.com

Technical details as well as toerance information can be requested seperately.